

# **Vaccine Manufacturing Using iCELLis<sup>®</sup> Reactor Systems Is Orchestrated by Lucullus<sup>®</sup> in a Validated Environment at IDT Biologika**

Pascal Vonlanthen, Bioprocess Specialist, Securecell AG  
Navjyot Waghmare, Automation Expert, IDT Biologika

## Vaccine Manufacturing Using iCELLis® Reactor Systems Is Orchestrated by Lucillus® in a Validated Environment at IDT Biologika

Written by:

Pascal Vonlanthen, Bioprocess Specialist, Securecell AG  
Navjyot Waghmare, Automation Expert, IDT Biologika

The demand for vaccines, cell and gene therapeutics, and biologic manufacturing in general has increased continuously over the past decade. In light of the COVID-19 pandemic, demand has reached new heights, with the requirement to produce billions of doses of COVID-19 vaccines and therapeutics in a short time. As an established CDMO for vaccines, gene and immune therapeutics, and other biologics, IDT Biologika has accumulated a wealth of experience in the manufacturing of such products – from process development and clinical manufacturing to commercial production.

For tomorrow’s technological progress, IDT employs innovative and proven technologies provided by leading software and hardware vendors. As part of our digitization strategy, we implemented Lucillus® software together with our partner Securecell.


The iCELLis® reactor systems (Pall) are an example of such innovative technology. These fixed bed bioreactors were designed for adherent cell culture applications such as viral vaccine and gene therapy vector production. In iCELLis® bioreactors, cells attach to a porous matrix which creates a huge surface area to optimize cultivation conditions and thus yield compared to commonly used hyper stack cell factories or microcarrier adherent cultures. Moreover, the iCELLis® technology is flexible from small scale to commercial volume, supporting an easier process scale-up.

Despite all the advantages the iCELLis® reactor systems have, they still operate as individual units and overarching process monitoring and control of multiple iCELLis® reactors and other peripheral devices is not standardly available. However, such overarching process monitoring and control capabilities were strongly desired by operators at IDT. A major complication to achieving overarching operation of the various scales (iCELLis® Nano, iCELLis® 500) is the different types of control units that control the vessels. The iCELLis® Nanos are controlled by Applikon my-controls, whereas the iCELLis® 500 use a very different automation




Exemplary image of an iCELLis 500 reactor empowered with Lucillus® at IDT in Braunschweig.

platform by Pall. Thus, standardized and centralized data logging, process control, and reporting were inherently difficult and could not be solved by either of the two hardware vendors. Furthermore, the operators at IDT had the desire to flexibly design their process control recipes (step-chains) themselves without dependency on the hardware manufacturer, which was often accompanied by long lead times and unnecessary financial efforts. In conclusion, Securecell AG was onboarded by Applikon BV for this integration project.



Example of an interactive piping and instrumentation diagram designed by IDT scientists to monitor and control the most crucial process parameters of their iCELLis systems®.



Network diagram of the client-server architecture across four labs at IDT Dessau encompassing 4 iCELLis® Nanos, 1 iCELLis® 500 and 4 iCELLis® 500+ reactors. Additionally, a pH sensor was integrated over a customized I/O box converting the analog signal to ethernet.

Navyjot Waghmare (*Automation Expert at IDT*) has led the subsequent integration project that resulted in Lucillus® empowered iCELLis® reactors in validated environments at IDT at its sites in Dessau and in Magdeburg supporting vaccine production to date. But first things first: The initial installation encompassed the integration of 4 iCELLis® Nanos and an iCELLis® 500 reactor in Lucillus® at IDT Dessau in 2017. Specifically engineered product improvements were also part of that implementation and exemplify Securecell's service promise and competence. One new feature that was implemented upon feedback from IDT was a configurable user interface allowing interactive piping and instrumentation diagrams (P&ID) to operate the equipment. This is a nice example of how proactive customers contribute to the targeted development of Lucillus® based on real use cases (from which other customers still profit today). In the following, Lucillus® was tested thoroughly and as a result has been operating in that validated production environment thus far.

In the subsequent years, Lucillus® was also adopted at the IDT development site in Magdeburg on another set of small-

scale iCELLis® Nanos. In addition, the initial installation in Dessau was extended by the integration of four additional iCELLis 500+, the newer generation of large-scale iCELLis® reactors, in September and November 2021, and now spans four labs in a robust, distributed client-server architecture.

With Lucillus®, IDT process engineers can now flexibly implement automated process control recipes to allow for quicker project adaptations. The operators can rely on unified monitoring and control interfaces built according to their needs. All data, including contextual information about the processes, is securely stored in a centralized database in a harmonized manner. Together, IDT and Securecell realized a performant bioprocess infrastructure with Lucillus® as the digital backbone, ensuring data integrity, greater efficiency, and increased capacity for IDTs offerings as a CDMO. This unique flagship implementation on fixed-bed production vessels will help IDT to manufacture and deliver high-quality and safe vaccines.

» *Every customer application is different; Lucillus® is extremely flexible and even when certain use cases are not yet covered, such as at IDT, product improvements will be evaluated by a customer-focused, solution-oriented team. It is these common journeys that enabled Lucillus® to grow into the powerful tool it is today.*

**Kay Eilers; Head of Business Development at Securecell AG**



»» We use the iCELLis® systems for many different processes and customer applications. Therefore, it was important to have a software solution like Lucullus®, which is flexible and can be easily adapted in a short time for new processes. The P&ID function implemented helps us to adapt the operator screen without dependence on system suppliers and to oversee the critical process parameters. The operators can rely on unified monitoring and control interfaces built according to their needs. All data, including contextual information about the processes, is securely stored in a centralized database in a harmonized manner.«

**Navjyot Waghmare, Automation Expert, IDT Biologika**



Both companies are now looking into handling Lucullus® life-cycle management of the installation in the validated space, which means going through an update process. Within the framework of the renewal, IDT will profit from the new, ever-improving, and increasing Lucullus® functionality. Specific examples are the REST API interface to facilitate data export to 3rd party data analysis and modeling tools, manual user interactions along the automated process or improved reporting to name just a few.

Within the last five years, we have integrated more than ten bioreactors at IDT in Lucullus® for GMP processes. In the future, the goal of both companies is to upgrade the Lucullus® software to the expert version, which will enable complete process digitalization from the planning until the process execution phase and further automate and digitalize GMP processes at IDT.

## References

1. <https://www.ema.europa.eu/en/medicines/human/EPAR/covid-19-vaccine-inactivated-adjuvanted-valneva>.
2. <https://valneva.com/press-release/valneva-and-idt-biologika-announce-collaboration-for-production-of-inactivated-covid-19-vaccine-vla2001/?lang=de>.

## About the company

In 1921, IDT Biologika started as a state institute with three employees. Today, 100 years later, IDT Biologika is a global biopharmaceutical contract development and manufacturing organization that specializes in the production of innovative live viral vaccines, viral vectors for gene & immune therapeutics, oncolytic viruses, virus-like particles and other sterile liquid or lyophilized biologics to improve human health worldwide. Our Analytical Development Department uses NGS for detecting adventitious viruses in the manufacturing processes. We validated our NGS analytics platform, designed and developed a bioinformatics workflow and validated it in conformance with GAMP 5 after qualifying the required IT infrastructure.